

MICROCOPY RESOLUTION TEST CHART NATIONAL BUREAU OF STANDARDS-1963-A

FTD-ID(RS)T-0803-82

FOREIGN TECHNOLOGY DIVISION



EXPERIMENTAL STUDIES OF A INVERTED CAVITY DYE LASER

by:

Wang Lianzhi, Yao Xinzo, Qiu Yuanwu



JUL 1 1 1984

A

Approved for public release; distribution unlimited.

84 07 10 119

EDITED TRANSLATION

FTD-ID(RS)T-0803-82

12 August 1982

MICROFICHE NR: FTD-82-C-001089

EXPERIMENTAL STUDIES OF A INVERTED CAVITY DYE LASER

By: Wang Lianzhi, Yao Xinzo, Qiu Yuanwu

English pages: 7

Source: Wuli, Vol. 10, Nr. 9, September 1981,

pp. 550-552

Country of origin: China Translated by: SCITRAN

F33657-81-D-0263

Requester: FTD/TQTD

Approve a for public release; distribution unlimited.

TMIS TRANSLATION IS A RENDITION OF THE ORIGINAL FOREIGN TEXT WITHOUT ANY ANALYTICAL OR EDITORIAL COMMENT. STATEMENTS OR THEORIES ADVOCATED OR IMPLIED ARE THOSE OF THE SOURCE AND DO MOT NECESSARILY REFLECT THE POSITION OR OPINION OF THE FOREIGN TECHNOLOGY DIVISION.

PREPARED BY:

TRANSLATION DIVISION FOREIGN TECHNOLOGY DIVISION WP-AFB, OHIO.

FTD -ID(RS)T-0803-82

Date 12 Aug 19 82

GRAPHICS DISCLAIMER

All figures, graphics, tables, equations, etc. merged into this translation were extracted from the best quality copy available.

Accession For

NTIS GRA&I
PLIC TAB
Unerpounced
Justification



A-1

EXPERIMENTAL STUDIES OF A INVERTED CAVITY DYE LASER Wang Lianzhi Yao Xinzo Qiu Yuanwu (Piysics Institute, Academia Sinica) Received August 26, 1980

The dye laser has been widely applied in many fields. The output /55(power of the ordinary pulsed optical pumping dye laser is of the order of several tens to several hundred watts, and its pulse width is of the order of microseconds. However, in many research areas, one would like to have dye laser output with high pulse power and narrow pulse width. Since the life time of the upper energy states is very short (e.g. the singlet upper energy level life time of Rhodamine 6G is 5.5x10⁻⁹ sec), it is very difficult to use the method of Q switch in the laser cavity to obtain a high powered pulse. On the other hand, an effective way to achieve high power and short pulse is the cavity inversion method. In this paper we introduced the experimental studies of the cavity inversion method used on the dye laser. The dye used is Rhodamine 6G; concentration is 10-4 mole. The inverted cavity pulse obtained has a wave length of 5900 Å and a pulse width of 11 millimicrosec. The output energy is 50 millijoules. The experimental results are basically in agreement with the theory. In this paper, we also analysed the key questions of getting cavity inversion laser output with high efficiency.

1. Basic Principles

The cavity inversion method has been used to obtain laser output with high efficiency and short pulse with results on several types of lasers. Preliminary investigation has also been carried out on the dye lasers.

The cavity inversion principle of the dye laser is shown in Figure 1. It is different from the ordinary laser resonance cavity. The two reflecting mirrors that form the resonance cavity are both total reflecting mirrors. Cavity inversion elements consisting of a




Figure 1.

Key: 1--dye medium; 2-- Pockels cell; 3-- polarization prism; 4--total reflecting mirror; 5--inverted cavity output

Pockels cell and a polarizing prism are installed in the cavity. Photons from the excited dye medium oscillate back and forth in the totally reflecting cavity. The optical pumping process continuously accumulates photons with hardly any laser output from the cavity. In this process, the Pockels cell and the polarizing prism guarantee that the resonance cavity is kept in a working state of high Q value. When the value of photon density in the cavity reaches a maximum, a half wave voltage is applied to the Pockels cell. As the polarized light passes through it, the polarization surface rotates by 90°, the photons are output in the direction of the perpendicular optical axis by the polarizing prism until all the photons in the cavity are exhausted. Hence, the output laser pulse width is

$$T = \frac{2L}{c} , \qquad (1)$$

where L is the cavity length and c is the speed of light. Before the potential on the Pockels cell is turned off, there will be no laser output once the photons in the cavity are exhausted because the resonance cavity is in a high output loss state. The dye is a 4 energy level system. To analyse the cavity inversion effect, -R.G. Morton [1], et al used a set of much simplified rate equations and obtained the approximate analytical expression for the photon number density in the cavity as

$$\theta = \frac{E^{D}}{E^{H}} = \frac{(2L/c)\phi_{p}^{D}}{\tau_{L}\phi_{p}^{D}}, \tag{2}$$

where E^{D} , E^{N} are the laser output energy (D represents the cavity inversion method and N the ordinary method).

photon number densities in the cavity, r_i is the laser pulse width for the ordinary method. Hence, other than suitably lengthening the cavity length L, increasing the cavity photon number density . in the cavity inversion method is a more important key in obtaining inverted cavity output of high efficiency. In the cavity inversion dye laser, the factors affecting the photon number density in the cavity are (1) the singlet state self-absorption of the dye medium. (2) triplet state laser absorption of the dye medium, (3) loss due /551 to the insertion of optical elements in the cavity such as the Pockels cell, the polarizing prism, etc. With the exception of loss factors (1) and (2) which are inherent to loss in the dye medium, the key question . of improving the inverted cavity efficiency is to reduce as much as possible the insertion loss of the laser. Reference [1] estimated that when the insertion loss was 0.1 (i.e., the percentage loss for hotons to travel once back and forth in the cavity), the photon number density in the cavity inversion method is dozens of times larger than that in the ordinary method. Thus, by carefully considering the insertion loss, it is entirely possible to obtain inverted cavity laser output with high efficiency.

2. Experimental Method and Results

In order to carry out the measurement of resonant fluorescence spectrum in the plasma diagnostic technique (such as sputtering impurity diagnosis) and to use a color center laser as optical pumping laser light source, we have investigated the cavity inversion laser with high power and short pulse. The experimental set up is shown in Figure 2.




Figure 2. Experimental set-up.

Key: 1--dye medium; 2--KDP crystal; 3--OE total output prism;
4--R2 (total reflecting); 5,7--diode; 6,8--oscilloscope; 9--power
meter; 10--R1=5m (total reflecting).

To improve high power laser output, it is very necessary to use short pulsed optical pumping light sources. The optical source we have chosen is a co-axial cylindrical lamp. The co-axial lamp tube consists of 2 quartz tubes. Ø of the inner tube are 15 and 13 mm, Ø of the outer tube are 19 and 17 mm. The tubes and the electrodes are co-axial in structure with a distance of 1 mm between the tubes. The discharge length is 35cm. The silvered focussing mirror is tubular with diameter 35 mm and length 300 mm. The concentration of Rhodamine 6G Ethyl alcohol dye solution is 10⁻⁴ mole. The 1 mm gap between the tubes is the discharge chamber filled with flowing pure Argon gas with pressure of 50-80 torr. Flowing distilled water passes between the focussing mirror and the outer tube to cool the co-axial tubes on one hand and to strengthen the ability of the co-axial lamp support on the other to receive strong current discharge without shattering. Since the dye container is long, both ends of the dye container are sealed into Brewster angle windows to prevent the production of self-oscillating laser light and to generate polarized laser light.

The resonance cavity of the laser i formed by 2 totally reflecting mirrors coated with multi-layers of dielectric films. The output side mirror is plane and the other is a concave spherical surface of radius of curvature 5m. This insures that along the long cavity, we have a structure of close semi-co-focussing. The Gauss light beam in the cavity is focussed in the vicinity of the polarizing prism to prevent the photons in the cavity leaking out. The Pockels cell is an electro-optic switch formed by KDP crystals. Under our conditions, the half wave potential is 9 K_V. The polarizing prism is a Glan OE total output prism with a hole diameter of 18 mm. This will reduce the insertion loss of the polarizing prism as much as possible. By experimental measurement, the transmissivity for the prism used is 88%.

The energy of the O light laser beam from the polarizing prism is recorded when it passes through a power meter. The pulse shape is displaye on an OK-19 model oscilloscope with a high current diode. We also recorded the optical pulse signal inside the cavity

from the small amount of photons leaking out through mirror R_2 . The entire discharge process block diagram is shown in Figure 3.




Figure 3. Synchronization block diagram of the discharge process.

Key: 1--gas-filled spark gap; 2--thyratron; 3--delayer; 4--hydrogen thyratron; 5--oscilloscope; 6--KDP crystal; 7--induction coil; 8--co-axial lamp.

We have obtained some results in our experiment. Figure 4 is the synchronous signal photographed from the Model OK-19 oscilloscope, indicating that the optical pulse signal in the cavity is well synchronized with the Pockels cell switch signal. Figure 5 gives the optical pulse signal in the resonance cavity. The pulse half width is about 0.9 µ sec. Figure 6 gives the inverted cavity laser output pulse signal obtained when the photon number censity in the cavity reaches peak value together with the synchronous Pockels cell potential. The laser wave length is 5900 Å, output is 50 millijoules, pulse width, 11 millimicrosecond. From the theoretical formula (1), _calculated with cavity length L=180cm, we obtained T=12 milli-microsecond. The experimental result is basically in agreement with theoretical calculations. From this calculation, the output laser power is found to be 4 megawatt. In the experiment, we also used resonance cavity with R_=00% reflectivity and measured the output signal of an /552 ordinary dye laser. Experiment indicates that the shape of the inverted cavity laser output is closely related to the correct switching time of the box.

Figure 6(a) and (c) show the inverted cavity output signal for leading and delaying switching pulse respectively. When the switching pulse precedes the maximal photon number density in the cavity, a definite multiple-peak phenomenor appears. When the switching pulse lags behind the maximal photon number density in the cavity, the inverted cavity laser pulse is significantly broadened. Only when the switching pulse is synchronized with the maximal photon number density in the cavity are we able to obtain a good output signal as shown in Figure 6 'b).

Summarizing the above statements, we see that by applying the cavity inversion method to dye lasers, it is possible to obtain high power and short pulsed laser output. The pulse width is completely determined by the cavity length. It may vary from a few milli-second to 100 milli-microsecond. The efficiency of the cavity inversion laser output is basically determined by the insertion loss of optical elements in the cavity. It is very important to construct low loss Pockels cell and polarizing prism.

It must be specifically pointed out that the high voltage open time of the Pockels cell must be much smaller than the inverted cavity pulse width. When the high voltage open time (i.e., the use time) is comparable to the inverted cavity pulse width, the cavity efficiency will be severely affected. We have developed a triode gas-filled co-axial spark gap composite switch as the high voltage pulse switch. The rising, leading edge of its pulse may reach 5 milli-microsecond or even as short as 1-2 milli-microsecond. This greatly increases the inverted cavity efficiency. In addition, if the dye laser is required to have tuned harmonic output, such harmonic tuning elements as the diffraction grating, Faraday rotator, standardizer, etc. may be applied to the inverted cavity dye laser.




Figure 4. A synchronous discharge signal. In this photograph the upper line is the optical pulse signal in the cavity of the inverted cavity dye laser; the lower line is the signal of the Pockels cell potential after voltage division.




Figure 5. Optical pulse signal in the resonance cavity 0.1 loss




Figure 6. a) Maximum photon number in the cavity before the switching pulse.

- b) Maximum photon number in the cavity synchronous to the switching pulse.
- c) Maximum photon number density in the cavity lagging behind the switching pulse.

Reference

[1] R.G. Morton et al., Appl. Opt., 17-20 (1978), 3268.

DISTRIBUTION CASE

DISTRIBUTION OFFICE OF STATEMENT

A205 LMWFN: A210 LMWW: B344 DIA/RTS-2C C043 USANITA C500 TRADOC C509 RALLISTIC RES LAB C510 RET LABS/AVRADXOM C513 ARRADXOM 1 C513 ARRADXOM 1 C535 AVRADCOM/TSARCOM 1 C539 TRASANA 1 C591 ESTC C619 MIA REDSTONE 1 D008 NISC 1 E053 HQ USAF/INET 1 E403 AFSC/INA 1 E404 AEDC/XOF 1 E408 AFWL E410 AD/IND 1 E429 SD/IND 1 P005 DOE/ISA/DDI P050 CIA/OCR/ADD/SD AFIT/LDE FTD CCN NIA/PIS NITS 1 NIA/PIS NITS 1 NISA/INST-44 NSA/INST-44 1 1 NSA/INS	OB WHINKLION	MICHOFICHE
B344 DIA/RTS-2C 9 C043 USAMITA 1 C500 TRADOC 1 C509 RALLISTIC RES LAB 1 C510 RET LABS/AVRADAOM 1 C513 ARRADAOM 1 C513 ARRADAOM 1 C535 AVRADAOM/TSARCOM 1 C539 TRASANA 1 C539 TRASANA 1 C691 ESTC 4 C619 MIA REDSTUNE 1 D008 NISC 1 E053 HQ USAF/INET 1 E403 AFSC/INA 1 E404 AEDC/DOF 1 E408 AFWL 1 E410 AD/IND 1 E429 SD/IND 1 E429 SD/IND 1 E429 SD/IND 1 FTD CCN 1 CCN CCN 1 CCN CCN 1 CCN CCN 1 CCN CCN CCN 1 CCN C	A205 DMARING	1
C043 USAMITA C500 TRADOC 1 C509 RALLISTIC RFS LAB C510 R&T LABS/AVRADAM 1 C513 ARRADYM 1 C535 AVRADCOM/TSARCOM 1 C539 TRASANA 1 C591 FSTC 4 C619 MIA REDSTONE 1 D008 NISC 1 E053 HQ USAF/INET 1 E403 AFSC/INA 1 E404 AEDC/DOF 1 E408 AFWL E410 AD/IND 1 E429 SD/IND P005 DOE/ISA/DDI P050 CIA/OUR/ADD/SD AFIT/LDE FTD CUN NIAS/PIS NIIS 2 ILINL/Cole L-389 NASA/NST-44	A210 1344VK	1
C500 TRADOC C509 BALLISTIC RES LAB C510 R6T LABS/AVRADXOM 1 C513 ARRADXOM 1 C513 ARRADXOM 1 C535 AVRADCOM/TSARCOM 1 C539 TRASANA 1 C591 FSTC 4 C619 MIA REDSTUNE 1 E008 NISC 1 E008 NISC 1 E403 AFSC/INA 1 E403 AFSC/INA 1 E404 AEDC/DOF 1 E408 AFWL 1 E410 AD/IND 1 E429 SD/IND 1 PO05 DOE/ISA/DDI PO05 CIA/OUR/ADD/SD AFIT/LDE FTD CUN NIA/PHS NIIS 1 NIIS 1 INII/Code L-389 NASA/NST-44	B344 DIA/RIS-2C	9
C509 BALLISTIC RES LAB 1 C510 R&T LABS/AVRADUOM 1 C513 ARRADUOM 1 C513 ARRADUOM 1 C513 ARRADUOM 1 C535 AVRADUOM/TSARUOM 1 C539 TRASANA 1 C591 FSTC 4 C619 MIA REDSTURE 1 D008 NISC 1 E053 HQ USAF/INET 1 E403 AFSC/INA 1 E404 AFDC/DOF 1 E408 AFWL 1 E410 AD/IND 1 E429 SD/IND 1 E429 SD/IND 1 E429 SD/IND 1 P005 DOF/ISA/DDI 1 P050 CIA/OUR/ADD/SD 2 AFIT/LDE 1 FTD CUN NIA/PRIS 1 NIIS 1 NIIS 1 NIIS 1 NIIS 2 IINL/Cole L-389 NASA/RST-44 1 C535 AVRADUOM 1 C515 AVRADUOM C515 AVR	CO43 USAMIIA	1
C510 R&T LABS/AVRADXOM C513 ARRADXOM C535 AVRADXOM/TSARXOM C539 TRASANA C591 FSTC C619 MIA REDSTONE D008 NISC E053 HQ USAF/INET E403 AFSC/INA 1 E404 AEDC/DOF E408 AFWL E410 AD/IND 1 E429 SD/IND 1 P050 CIA/OUR/ADD/SD AFIT/LDE FTD CCN NIA/PIS NIIS 1 INIA/PIS INIIS 1 INIA/Code L-389 NASA/NST-44	C500 TRADOC	1
C510 R&T LABS/AVRADXOM C513 ARRADXOM C535 AVRADXOM/TSARXOM C539 TRASANA C591 FSTC C619 MIA REDSTONE D008 NISC E053 HQ USAF/INET E403 AFSC/INA 1 E404 AEDC/DOF E408 AFWL E410 AD/IND 1 E429 SD/IND 1 P050 CIA/OUR/ADD/SD AFIT/LDE FTD CCN NIA/PIS NIIS 1 INIA/PIS INIIS 1 INIA/Code L-389 NASA/NST-44	C509 BALLISTIC RFS LAB	1
C513 ARRADYOM C535 AVRADCOM/TSARCOM C539 TRASANA C591 FSTC 4 C619 MIA REDSTONE D008 NISC I E053 HQ USAF/INET E403 AFSC/INA 1 E404 AEDC/LOF E408 AFWL 1 E410 AD/IND 1 E429 SD/IND 1 P050 DDE/ISA/DDI P050 CIA/OUR/ADD/SD AFIT/LDE FTD CCN NIA/PHS NIIS 1 NIIS 1 INIA/PHS NIIS 1 INIA/Code L-389 NASA/NST-44		• 1
C539 TRASANA C591 FSTC 4 C619 MIA REDSTONE D008 NISC I E053 HQ USAF/INET E403 AFSC/INA E404 AEDC/DOF E408 AFWL E410 AD/IND E429 SD/IND P005 DOE/ISA/DDI P050 CIA/OUR/ADD/SD AFIT/LDE FTD CUN NIA/PHS NIIS S IINL/Code L-389 NASA/NST-44		1
C591 FSTC C619 MIA REDSTUNE D008 NISC I E053 HQ USAF/INET E403 AFSC/INA E404 AEDC/DOF E408 AFWL I E410 AD/IND E429 SD/IND P005 DOE/ISA/DDI P050 CIA/OCR/ADD/SD AFIT/LDE FTD CUN NIA/PIS NIIS I INIA/PIS NIIS I INIA/CO/Je L-389 NASA/NST-44	C535 AVRADCOM/TSARCOM	1
C619 MIA REDSTUNE D008 NISC E053 HQ USAF/INET E403 AFSC/INA E404 AEDC/DOF E408 AFWL E410 AD/IND E429 SD/IND P005 DOE/ISA/DDI P050 CIA/OCR/ADD/SD AFIT/LDE FTD CUN NIA/PIS NIIS 1 NASA/NST-44	C539 TRASANA	1
D008 NISC E053 HQ USAF/INET E403 AFSC/INA E404 AEDC/DOF E408 AFWL E410 AD/IND E429 SD/IND P005 DOE/ISA/DDI P050 CIA/OCR/ADD/SD AFIT/LDE FTD CUN NIA/PIS NIIS 1 NASA/NST-44	C591 FSTC	· 4
E053 HQ USAF/INET E403 AFSC/INA E404 AEDC/DOF E408 AFWL E410 AD/IND E429 SD/IND E429 SD/IND P005 DOE/ISA/DDI P050 CIA/OCR/ADD/SD AFIT/LDE FTD CUN NIA/PIS NIIS 1 NIA/PIS 1 NASA/NST-44	C619 MIA REDSTONE	· 1
E403 AFSC/INA E404 AEDC/DOF E408 AFWL E410 AD/IND E429 SD/IND P005 DOE/ISA/DDI P050 CIA/OUR/ADD/SD AFIT/LDE FTD CUN NIA/PHS NIIS 1 NASA/NST-44	D008 NISC	1
### F404 AFDC/IXOF 1	E053 HQ USAF/INET	1
E408 AFWL 1 E410 AD/IND 1 E429 SD/IND 1 P005 DOE/ISA/DDI 1 P050 CIA/OUR/ADD/SD 2 AFIT/LDE 1 FTD 1 CUN 1 NIA/PIS 1 NIIS 2 IINL/Code L-389 1 NASA/NST-44 1		· 1
E410 AD/IND 1 E429 SD/IND 1 P005 DOE/ISA/DDI 1 P050 CIA/OUR/ADD/SD 2 AFIT/LDE 1 FTD 1 CUN 1 NIA/PIS 1 NIIS 2 IINL/Code L-389 1 NASA/NST-44 1	F404 AFDC/DOF	1
E429 SD/IND 1 P005 DOE/ISA/DDI 1 P050 CIA/OUR/ADD/SD 2 AFIT/LDE 1 FTD 1 CUN 1 NIA/PIS 1 NIIS 2 IINL/Code L-389 1 NASA/NST-44 1	E408 AFWL	1
P005 DOF/ISA/DDI 1 P050 CIA/OUR/ADD/SD 2 AFIT/LDF: 1 FTD 1 CUN 1 NIA/PIS 1 NIIS 2 IINL/Cole L-389 1 NASA/NST-44 1	E410 AD/IND	1
P050 CIN/OUR/ADD/SD 2 AFIT/LDE 1 FTD 1 CUN 1 NIA/PHS 1 NIIS 2 IINL/Cole L-389 1 NASA/NST-44 1	E429 SD/IND	1
AFIT/LDE 1 FTD	P005 DOE/ISA/DDI	1
AFIT/LDE 1 FTD	POSO CIA/OUR/ADD/SD	2
CCN 1 NIA/PHS 1 NIIS 2 IINL/Cole L-389 1 NASA/NST-44 1		1
NIA/PIIS 1 NIIS 2 IINL/Cale L-389 1 NASA/NST-44 1	FTD	
NIIS 2 IINL/Cole L-389 1 NASA/NST-44 1	CON '	. 1
11NL/Cale L-389 1 NASA/NST-44 1	NIA/PIIS	1
NASA/NST-44		2
	IINL/Cale L-389	$\overline{1}$
	NASA/NST-44	ī
		. 2

